

## Available formats:

Container 60 ml

Container 150 ml

Container with 500 ml gun



**HYPOCLORINE**<sup>®</sup>  
Ácido Hipocloroso



*Healthy skin*

# HYPOLCLORINE

## SKIN CARE

### COMPOSITION

**Hypocloline skin care HYDROGEL:** Ionized Water, Hypochlorous Acid, Sodium Chloride, Sodium Fluorosilicate Gel, Magnesium, Monosodium Phosphate.

**Hypocloline skin care LIQUID:** Ionized Water, Hypochlorous Acid, Sodium Chloride.

Veterinary Formula intended for cleaning and hygiene of the skin of animals. It has a broad spectrum antimicrobial activity, destroying 99.9% of germs. Active against Gram (+), Gram (-) bacteria, viruses, fungi and yeasts. Eliminates odors by acting on the germs that generate them. neutral pH.

How to use **HYPOLCLORINE SKIN CARE** Veterinary formula for skin:

- Remove excess hair from the area. Cut the hair if necessary.
- Apply HYPOLCLORINE SKIN CARE generously to the area, according to your veterinarian's advice.
- In the area to be treated, cover it completely with **HYPOLCLORINE SKIN CARE**.
- When applying **HYPOLCLORINE SKIN CARE**, do not use any other topical product at the same time.
- If you consider it necessary, cover the area with a bandage or dressing, apply HYPOLCLORINE SKIN CARE generously on the bandage or dressing. No need to rinse or dilute.

### TARGET SPECIES:

Use in animals in general

### STORAGE:

Store at room temperature out of direct sunlight and heat. Do not allow to freeze. After use, close the container.