

Available formats:

Container 60ml



HYPOCLORINE[®]
Ácido Hipocloroso



Healthy ear

HYPOCLORINE

EYE CARE

COMPOSITION

Hypochlorine ear care HYDROGEL: Ionized Water, Hypochlorous Acid, Sodium Chloride, Sodium Fluorosilicate Gel, Magnesium, Monosodium Phosphate.

Hypochlorine ear care LIQUID: Ionized Water, Hypochlorous Acid, Sodium Chloride.

Veterinary Formula intended for cleaning and hygiene of the skin of animals. It has a broad spectrum antimicrobial activity, destroying 99.9% of germs. Active against Gram (+), Gram (-) bacteria, viruses, fungi and yeasts. Eliminates odors by acting on the germs that generate them. neutral pH.

How to use **HYPOCLORINE EAR CARE** Veterinary formula for external ear:

- Apply generously in the outer ear.
- The product is ready to be used without diluting or rinsing.

TARGET SPECIES:

Use in animals in general

STORAGE:

Store at room temperature out of direct sunlight and heat. Do not allow to freeze. After use, close the container.